

PROXY

Automatic Endoscope Reprocessing (AER) Chemistry System*

Key Benefits:

- ✓ Australian made
- ✓ Extremely cost efficient

Insist on Whiteley when ordering your Endoscopy products



PROXY A

Anticorrosive Additive

An additive for use with PROXY P.
ARTG No: 274095

2 x 5 Litres

Product Code: 2x5L - 160622



PROXY C+

Concentrated Pre-disinfectant Cleaner

An acidic pre-disinfectant cleaner, compatible with PROXY automatic endoscope reprocessing chemistry system. ARTG No: 125527

2 x 5 Litres

Product Code: 2x5L - 160642



PROXY P

High Level Disinfectant & Sterilant

A 5% peracetic acid based high level disinfectant and sterilant.

For use with PROXY A.

ARTG No: 274095

2 x 5 Litres

Product Code: 2x5L - 160612



PROXY CLN

Multi-Enzyme Detergent

Multi-Enzyme, low foaming detergent compatible with PROXY automatic endoscope reprocessing chemistry system. ARTG No: 125529

2 x 5 Litres

Product Code: 2x5L - 160652



PROXY L

70% Ethanol Solution

70% solution of Ethanol in water.

ARTG No: 127945

4 x 2 Litres

Product Code: 4x2L - 160632

*Suitable for use in Soluscope 3, and Sprint PAA AERs, please contact your Whiteley representative prior to first use.

ABN No. 24 000 906 678
PO Box 1076, North Sydney NSW 2059
Product Support Hotline: 1800 833 566
www.whiteley.com.au

NZBN: 9429041224243
PO Box 22-519, Otahuhu, Auckland 1640
Product Support Hotline: 0800 257 352
www.whiteley.com.nz

Whiteley 

Protecting lives and livelihoods since 1933