

## Laboratory Spill Response Kit

SKU: Laboratory Spill Response Kit (ZTSSLNK)

Link: <https://www.spillstation.com.au/spill-kits-lab-med/laboratory-spill-response-kit/>



## Product Information

---

### Absorbing Capacity: 12 Litre

The ZTSSLNK laboratory spill neutralising kit is designed to neutralise and absorb a broad range of acids, alkalis or base, hydrocarbons and toxic substances found in the laboratory. These kits also contain basic PPE and spill disposal tools.

Not suitable for hydrofluoric acid.

## Product Description

---

The acid and alkali or base neutralising polymers are formulated to move Class 8 substances towards neutral and at the same time, convert the spill into an easily managed solid. The multi-purpose absorbent powder is suitable for Class 6, Class 3 and biohazard absorption.

This spill kit comes complete with warning signs, instructions, waste disposal bags and basic personal protection equipment but additional equipment such as extra absorbents or organic vapour respirators can be included if required.

ZTSSLNK Contents:

- Alkali Neutralising Absorbent
- Acid Neutralising Absorbent
- Multi-Purpose Absorbent

- MAX60 Absorbent Sheets
- Protective Disposable Apron
- Disposable Nitrile Gloves
- Overshoes and Hairnet
- Face Mask and Eye Shield
- Scoop and Scraper
- Clinical Waste Bags
- "Caution Spill" Sign
- Labelled 10L Pail
- Instructions

Additional Personal Protection Equipment is available if required.

[/vc\_column][/vc\_row]