



## MEDI-DRY CONCENTRATE

#### **RINSE AID**

AUSTRALIAN MEDICAL DEVICE CLASS 1
ARTG No. 125529

## **BENEFITS:**

- Ideal for use in hard or soft water
- Concentrated formula
- Class 1 Medical Device
- · Promotes rapid drying of Instruments and Utensils
- Ultra low foaming

MEDI-DRY CONCENTRATE is a liquid rinse aid designed for use as a water dispersant in automatic washer/disinfectors. The concentrated formula is suitable for use in hard and soft water.

The primary function of MEDI-DRY CONCENTRATE is to disperse the rinse water to ensure the rapid drying of instruments and utensils. MEDI-DRY CONCENTRATE uses readily biodegradable ingredients.

### **MAJOR ADVANTAGES:**

## **Guaranteed quality:**

Use of only highly efficient ingredients ensures consistent quality performance, free from rejects. MEDI-DRY CONCENTRATE produces fast spot and streak free drying of all washed articles.

### **Consistent destaining:**

Excellent rinsibility gives the brightest possible results with less water spotting and streaks after heavy loads.

## **Better machine performance:**

MEDI-DRY CONCENTRATE acts as a mild foam suppressant in all single tank machines.

#### **DIRECTIONS FOR USE:**

Dispense or supply the washer/disinfector at a rate of 0.2 – 0.4mL per Litre of water.

Medi-Dry Concentrate - December 2011 1





# **TECHNICAL BULLETIN**

#### **CHEMICAL COMPOSITION AND PROPERTIES:**

Surfactant System:.... Non-ionic Colour (Concentrated):..... Blue Liquid Colour (Dilute solutions):.... Colourless Odour:.... **Odourless** pH (Concentrate):.... 6.5-8.5 pH (1% solution):.... 7.0-8.0 **Completely Soluble** Solubility:.... Rinsibility:.... Excellent, residue free Specific Gravity: 1.016 Foam Properties:.... Low foaming Dangerous Goods Rating:..... None Shelf Life: 2 years

**SAFETY DATA:** Information is listed on the Material Safety Data Sheet.

**PACK SIZE:** 2 x 5 Litre Wide Neck

#### **ENVIRONMENTAL STATEMENT:**

Whiteley Corporation manufactures products that comply with Australian and International environmental policy. Whiteley Corporation aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley Corporation products are supplied in 100% recyclable containers. Whiteley Corporation and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

© This document remains property of Whiteley Corporation Pty Ltd.

Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.

Medi-Dry Concentrate - December 2011 2

