

Solutions for Sterile Processing



Insist on Whiteley Medical when ordering your Sterilising products

Providing quality patient care can only be successfully achieved if the surgical instruments and equipment used to treat them are clean.

MACHINE DETERGENT

Whiteley Medical's washer disinfectant detergents are suitable for any brand of machine in the CSSD. Formulated using cutting edge technology, Whiteley Medical's machine detergents produce superior cleaning results on medical instruments.



Pack Size: 5L

GENESIS

Multi Function Machine Detergent

- Concentrated formula with minimal foaming.
- Cost savings with lower dosage and less inventory required.
- Compatible with commercially available simulated soil tests (including Browne's, TOSI and Edinburgh Soil Tests).
- Complies to requirements of ISO 15883-1:2006.
- Excellent materials compatibility with medical instruments.
- Suitable for use in hard and soft water.
- Excellent rinsability.
- Contains corrosion inhibitor.
- ARTG No: 125529.

GENESIS is a pH neutral, multi-function, enzymatic machine detergent developed for the automated cleaning of medical instruments. Its revolutionary formula utilises complex technology to produce superior cleaning efficacy on medical instruments with minimal foaming.

MACHINE RINSE AIDS



Pack Size: 5L

MEDIRINSE

Rinse Aid

- Disperses water to ensure the rapid drying of instruments and utensils.
- In built scale inhibitor.
- Ideal for use in hard or soft water.
- Concentrated formula.
- ARTG No: 125529.



Pack Size: 5L, 15L

ARID

Mildly Acidic Rinse Aid

- Acidic buffering system – reduces the need for descaling
- Low foaming – free rinsing
- Contains a corrosion inhibitor
- Phosphate free
- Dual Sequesterant System: prevents metallic deposition and staining
- Low viscosity – suitable for use in most washer/ disinfectors
- ARTG No:125527.

DESCALERS & NEUTRALISERS



Pack Size: 5L

MEDAC

Acidic Detergent

- Brightens stainless steel and glass instruments and utensils.
- Can be used for surface brightening of stainless steel benches and trolleys.
- Removes scale and soil deposits in sterilisers, washer disinfectors, ultrasonic machines and anaesthetic apparatus decontaminators.
- Assists in the removal of surface rust.
- Can be used to brighten instruments in an ultrasonic machine.
- ARTG NO: 125527.



Pack Size: 5L, 15L

PASAC

Scale Remover – Neutraliser

- Specifically formulated to prevent scale.
- Ideal for use in hard or soft water.
- Phosphate free.
- Liquid neutralising agent based on organic acids, guaranteeing rinsing of all alkali residues.
- ARTG No: 125527.

Whiteley Medical has a complete range of innovative products required for central sterile service departments.

MANUAL & ULTRASONIC DETERGENTS

Surgical instruments are often very expensive and delicate. The chemicals used to clean these instruments need to be specifically formulated to effectively clean without harming the instruments or users.



Pack Size: 500mL, 5L

MEDIZYME

Proteolytic Enzyme Cleaner

- Immediately effective on organic soils including blood, body tissue and carbohydrates.
- Highly active multi-enzyme formulation ensures superior results.
- Ideal for manual or ultrasonic cleaning.
- Free rinsing.
- ARTG No: 125529.



Pack Size: 1L, 5L, 15L

MATRIX®

Biofilm Remover

- Matrix® is the world's leading biofilm remover!
- Matrix® has been demonstrated to remove conventional biological soils and bacterial biofilms.
- Cleaning of medical devices has been identified as the critical first step in any effective instrument reprocessing system.
- Matrix® is the first documented and proven solution to the complex problem of validated cleaning of modern medical devices.
- ARTG No: 125529.



Pack Size: 1L, 5L

CLINIMAX®

Instrument & Equipment Detergent

- Concentrated, non-hazardous formula as per Safework Australia criteria
- Odourless
- Low foaming to improve visibility during cleaning processes
- Free rinsing
- Effective in hard and soft water
- Contains a corrosion inhibitor to protect quality instruments
- ARTG No: 125528.



Pack Size: 5L, 15L

SONIDET®

Medical Equipment & Instrument Detergent

- pH neutral – safe for use on all medical devices.
- Specifically designed for both ultrasonic and manual cleaning operations.
- Free rinsing and bacteriostatic.
- Odourless, non-staining and is a safe detergent for the removal of hospital and surgical soils.
- ARTG No: 125529.

SPECIAL PURPOSE



Pack Size: 500mL, 5L

MR STEEL

Stainless Steel Polish

- Water-based emulsion for the cleaning and protection of stainless steel surfaces.
- Restores stainless steel to as-new lustre.
- ARTG No: 187159

Looking for our Safety Data Sheets?

Whiteley Corporation provides Safety Data Sheets (SDS) for all products. Search for SDS on our website at: www.whiteley.com.au/SDS





Insist on Whiteley Medical when ordering your Sterilising products

HEALTH STATEMENT

Whiteley Medical makes every effort to ensure that its products are always safe to use, however products can affect individuals differently and the company cannot predict individual responses. Therefore careful reading of the SDS is essential prior to first use.

ENVIRONMENTAL STATEMENT

Whiteley Medical manufactures products that comply with Australian and International environmental policies.

DISTRIBUTED BY:



ABN No. 24 000 906 678
PO Box 1076, North Sydney NSW 2059
Product Support Hotline: 1800 833 566
www.whiteley.com.au

NZBN: 9429041224243
3-5 Kahu Street, Otahuhu, 2024
Product support Hotline: 0800 257 352
www.whiteley.co.nz